


# HORIBA JOBIN YVON

## DSS - SERIES PHOTODIODE / RECEIVER



**Application Note** This unit is a high performance photodiode/receiver operated with a thermoelectric cooler for stabilization/cooling with a dual gain FET input transimpedance amplifier. The output voltage is proportional to the input signal current:  $V_{out} = I_{signal} \cdot R_f$ . The PD/AMP is a DC coupled dual gain system. Care should be taken in shielding the unit from stray light during operation to prevent saturation of the amplifier (and potential failure).


### SPECIFICATIONS

Detector Type	3 mm ex-InGaAs Photodiode	
Operating Temperature - °C	25 @ $I_{tec} = 0.0 A$	-30 @ $I_{tec} = 0.6 A$
Operating Wavelength - $\mu m$	1.0 - 2.1	1.0 - 2.1
Responsivity - V/W @ pk	$1.2 \times 10^7 / 10^6$	$1.2 \times 10^7 / 10^6$
Noise - V/Hz <sup>1/2</sup> @ 100 Hz	$1.8 \times 10^{-5} / 10^{-6}$	$4.0 \times 10^{-6} / 10^{-7}$
NEP - W/Hz <sup>1/2</sup> @ 1.95 $\mu m$	$< 1.5 \times 10^{-12}$	$< 3 \times 10^{-13}$
Bandwidth (-3dB) - Hz	DC - 2k	DC - 2k
Power Requirements	+/- 9 VDC to +/- 15 VDC	
Connections	BNC signal output. Shielded power cable terminated with a DB-9 connector directly couples the unit with the PS/TC -1 Low Noise Power Supply.	

## ***DSS-SERIES SPECTRAL RESPONSE***

### **DSS-IGA1.9A PHOTODIODES**


**SPECTRAL RESPONSE - IGA1.9 AMBIENT-series**



## ***DSS-SERIES SPECTRAL RESPONSE***

### **DSS-IGA1.9TC PHOTODIODES**

**SPECTRAL RESPONSE - IGA1.9 TC-series**





### OPERATING THE DSS-SERIES PHOTODIODE/RECEIVER

**POWER SUPPLY:** A bipolar power supply is required,  $\pm 6\text{VDC}$  to  $\pm 15\text{VDC}$ , 20mA. This means a +V, central/common ground and a -V connection - 3 wires total, to pins 6, 7, & 8 on the D-sub connector. The power supply pins should be bypassed physically close to the amplifier module. Double check wiring prior to turning on power. Improper /reverse wiring will damage the unit.

**GAIN SELECT:** The unit is supplied with a switch which provides a 10:1 HI/LO gain function. HI Gain is the up position on the switch; LO Gain is the down position. Consult the individual data sheet for specific values.

**AMBIENT LIGHT:** Because of the high gains involved, the unit must be shielded from ambient background light during operation. Measurement errors and/or saturation can result from improper shielding.

**OUTPUT CONNECTION:** The signal output is thru a BNC connector (or BNC terminated cable in the case of the 2-color units) located on the back of the module.